

**IPL Apollo System – Intense Pulsed Light system provides safe and effective solutions for hair removal, facial rejuvenation and problem skin problems.**

**Treatments of the following 3 categories;**

The APOLLO is a high powered intense pulsed light hair removal and photo rejuvenation system provides safe and effective treatments of the following 3 categories;

**1) Hair Removal to treat the following:**

Hair removal skin types I to VI

**2) Photo-Rejuvenation treatments to treat the following:**

Lines and wrinkles  
Age spots  
Sun spots and freckles  
Sun damaged skin  
Active acne clearance  
Pigmentation marks

**3) Vascular lesions to treat the following:**

Redness of the face & Red Neck  
Reduces and helps to clear rosacea.  
Decreases pore size  
Stimulates collagen & elastin



**TECHNICAL SPECIFICATIONS**

- **Intense Pulsed Light**
- **Light source: xenon lamp**
- **Energy per single shot: 80-150J**
- **Fluence max: 40J/cm<sup>2</sup>**
- **Pulse width: 5-130 ms**
- **Emission mode: single shot, burst - train 2-6 shots**
- **Standard filter: 530-1200 nm**
- **Cooling system: water (hand-piece included)**

Intense Pulsed Light (IPL) the energy for a single shot is between 80 and 150J, the pulse width is from 5 to 130 ms, this enables you to perform a diverse range of treatments in complete safety. IPL technology works particularly well on pale skins and larger areas.

**COMPLETE PACKAGE**

The Apollo is supplied as a complete package and consists of the following as standard:-

- IPL system in stand alone case
- Standard IPL hand-piece size 50 x 10 mm
- Water Cooled detachable hand-piece
- Filter wavelength 530-1200 nm
- Variable Pulse technology
- Touch screen and preset parameters for IPL Hair Removal & Photo Rejuvenation
- Manual parameter changes
- Goggles for eye protection
- 12 Month Warranty – direct from manufacturer
- Point of Sale Material
- Helpline for treatment queries
- Training

**Accredited Training is provided with the Apollo System  
Includes 12 months Insurance**